

Thermo Scientific explosion-proof refrigerators, freezers and combination units are designed to be used in areas where ignitable materials are normally present.

A range of models, particularly suited for safe storage in hazardous environments, are available with sizes from undercounter to full height freestanding.

Not available outside North America.

## Thermo Scientific Explosion-Proof Laboratory Refrigerator and Freezers



5.6 cu. ft. Explosion-Proof Refrigerator  
Model: 3557-DB

- Sparkfree interior and exteriors reduce the risk of explosion
- All models are manual defrost
- Tough white exterior
- CFC- and HCFC-free refrigeration system and insulation
- Designed for use in Class 1, Division 1, Group C and D environments
- Plug-free cords (units must be hardwired to wall)

**Warranty:** 13 months parts and labor

**Certifications:** CSA and or UL



5 cu. ft. Explosion-Proof Freezer  
Model: 3557-2



20.6 cu. ft. Explosion-Proof Freezer  
Model: 3552A



10 cu. ft. Explosion-Proof  
Combination Unit  
Model: 3551-DB

**Thermo Scientific Explosion-Proof Laboratory Refrigerators and Freezers** (Not Available outside North America)

**Refrigerators (Temperature Range 1° to 10°C)**

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves (interior)	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3557-DB	5.6	manual	120V / 60Hz	4 adjustable	yes	16 x 20 x 27.5**	29.5 x 23.5 x 33.5	110
3650A*	21	manual	120V / 60Hz	4 adjustable	yes	16.5 x 27 x 55	30.5 x 32 x 74.5	230
3566A	21	manual	120V / 60Hz	4 fixed	yes	16.5 x 27 x 55	31.75 x 32 x 74.5	230

\* Corrosion-resistant model, features stainless steel interior door liner and stainless steel shelves

\*\* 10" (for bottom 8"); 16" (for top 19.5")

**Freezers (Temperature Range -12° to -20°C)**

Model No.	Capacity	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3557-2	5	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 18 x 28.5	29.75 x 21 x 33.5	84
3552A	20.6	manual	115V / 60Hz	4 fixed	yes	28.5 x 32 x 70.8	26.5 x 32 x 73	200

**Combination Refrigerator/Freezers (Temperature Range 1° to 10°C Ref., -12° to -20°C Frz.)**

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3557-4	5.6	manual	120V / 60Hz	Ref: 3 fixed 3 door	no	18 x 15 x 21 (ref.) 11 x 15 x 5 (frz.)	25.5 x 21.25 x 32.75	78
3551-DB	10	manual	120V / 60Hz	Ref: 4 door, 3 adjustable; Frz: 1 fixed	no	13.5 x 19.5 x 38.5 (ref.) 15 x 17.5 x 11.5 (frz.)	31.5 x 23.5 x 60	150



**5.6 cu. ft. Explosion-Proof Combination Unit  
Model: 3557-4**



© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**North America:** USA/Canada +1 866 984 3766 (866-9-THERMO)

[www.thermoscientific.com/cold](http://www.thermoscientific.com/cold)

**Europe:** Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

**Asia:** China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,

Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940 or +33 2 2803 2180

DSCSEXPLPRF 0111