

Thermo Scientific flammable storage refrigerators, freezers and combination units are designed to be used for storage of volatiles where a spark-free internal cabinet environment is required.

A range of models particularly suited for storage of flammable materials with sizes from small undercounter to full height freestanding.

Not available outside North America.

## Thermo Scientific Flammable Storage Laboratory Refrigerator and Freezers



5.6 cu. ft. Flammable Storage Combination Unit  
Model: 3556-6

- Spark-free interiors to reduce risk of internal explosion
- All units are manual defrost
- Tough white exterior (except model no. 3575-10)
- CFC- and HCFC-free refrigeration system and insulation
- Designed to meet NFPA 45 and 99 along with OSHA article 1910.307

**Warranty:** 13 months parts and labor

**Certifications:** CSA and or UL



50 cu. ft. Flammable Storage Refrigerator  
Model: 3575-10



21 cu. ft. Flammable Storage  
Refrigerator  
Model: 3566-10A



10 cu. ft. Flammable  
Storage Combination Unit  
Model: 3551-10

**Thermo Scientific Flammable Storage Laboratory Refrigerators and Freezers** (Not Available outside North America)

**Refrigerators (Temperature Range 1° to 10°C)**

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves (interior)	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3556-DB	5.6	manual	120V / 60Hz	4 adjustable	yes	16 x 20 x 27.5**	28.75 x 23.5 x 33.5	110
3560-10A*	21	manual	120V / 60Hz	4 adjustable	yes	16.5 x 27 x 55	30.5 x 32 x 74.5	230
3566-10A	21	manual	120V / 60Hz	4 fixed	yes	16.5 x 27 x 55	31.75 x 32 x 74.5	230
3575-10†	50	manual	115V / 60Hz	8 adjustable	yes	26.5 x 45.75 x 59	34.25 x 52 x 77.25	400

\* Corrosion-resistant model, features stainless steel interior door liner and stainless steel shelves

\*\* 10" (for bottom 8"); 16" (for top 19.5")

† Features digital readout; stainless steel interior and exterior

**Freezers (Temperature Range -12° to -20°C)**

Model No.	Capacity	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3556-4	5	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 18 x 28.5	31.4 x 21 x 33.5	84
3552-10A	20.6	manual	115V / 60Hz	4 fixed	yes	28.5 x 32 x 70.10	26.5 x 32 x 73	200

**Combination Refrigerator/Freezers (Temperature Range 1° to 10°C Ref., -12° to -20°C Frz.)**

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3556-6	5.6	manual	120V / 60Hz	Ref: 3 fixed, 3 door	no	18 x 15 x 21 (ref.) 11 x 15 x 5 (frz.)	25.5 x 21.25 x 32.75	78
3551-10	10	manual	120V / 60Hz	Ref: 4 door, 3 adjustable Frz: 1 adjustable	no	13.5 x 19.5 x 38.5 (ref.) 15 x 17.5 x 11.5 (frz.)	30 x 23.5 x 60	150



**5 cu. ft. Flammable Storage Freezer**  
Model: 3556-4



**5.6 cu. ft. Flammable Storage Refrigerator**  
Model: 3556-DB



© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**North America:** USA/Canada +1 866 984 3766 (866-9-THERMO)

[www.thermoscientific.com/cold](http://www.thermoscientific.com/cold)

**Europe:** Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

**Asia:** China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,

Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940 or +33 2 2803 2180

DSCSFLAMM 0111