

Edition 04/2008
Replace Ed. 04/2005

Malt Extract R2 19595

Definition

Malt Extract R2 is manufactured from malt by a purification process.

Description

Fine cream powder easily soluble in water.

Malt Extract R2 contains a mix of carbohydrates (mainly maltose) and growth factors.

Use

Source of carbohydrates recommended in media for:

- Analytical microbiology
- Industrial fermentation.

Physico-chemical characteristics

	Standard
Solubility in water at 3 %	Complete
pH (3 % solution)	4.8 - 5.8
Loss on drying	≤ 6.0 %
Reducing sugars (as maltose)	≥ 60.0 %
Residue on ignition	≤ 4.5 %
Chloride (as NaCl)	≤ 1.0 %

Microbiology

	Standard
Total aerobic microbial count	≤ 10 000 /g
<i>Escherichia coli</i>	Absence / 10 g
<i>Salmonella</i>	Absence / 25 g
Yeasts and moulds	≤ 20 /g



Organotechnie® S.A.S.

27, avenue Jean Mermoz
93120 La Courneuve, France
Tél : +33 (0) 1 49 92 87 50
Fax : +33 (0) 1 49 92 87 51

e-mail : info@organotechnie.com
web : <http://www.organotechnie.com>



Les informations contenues dans ce document, données à titre indicatif, sont conformes à nos connaissances actuelles.

Il est de la responsabilité des utilisateurs de conduire leurs propres tests pour déterminer les conditions d'utilisations selon leurs usages spécifiques.

The information contained in this publication is based on our own research and development work and is to the best of our knowledge true and accurate.

Users should, however, conduct their own tests to determine the suitability of our products for their own specific purposes.

Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents.

Malt Extract R2 / 19595

Edition 04/2008
Replace Ed. 04/2005



Organotechnie S.A.S.

Documentation

The certificate of analysis and the sanitary certificate are supplied with each delivery.

Packing and storage

25 kg net corrugated board box with inner polyethylene bags.

Upon request: 5 kg plastic drum.

Keep in original packaging closed when not in use,
at room temperature in a dry area.

Hygroscopic product.

Best before: 5 years.

Health and safety information

Dusty powder.

Avoid inhalation.

Les informations contenues dans ce document, données à titre indicatif, sont conformes à nos connaissances actuelles.

Il est de la responsabilité des utilisateurs de conduire leurs propres tests pour déterminer les conditions d'utilisations selon leurs usages spécifiques.

The information contained in this publication is based on our own research and development work and is to the best of our knowledge true and accurate.

Users should, however, conduct their own tests to determine the suitability of our products for their own specific purposes.

Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents.